CARESITE® Smallbore Y-Extension Set

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en
CARESITE® Smallbore Y-Extension Set with two CARESITE Luer Access Devices and Spin-Lock® Connector
INTENDED PURPOSE: CARESITE Extension sets are sterile, single use devices for direct injection, intermittent infusion, continuous infusion or aspiration of fluids, medications, blood and blood products. The device is accessed via the CARESITE luer access device upon insertion of a male luer.
TARGET PATIENT GROUP: Extension sets can be used for all patients who require IV therapy.
INTENDED USER: Extension Sets are used by trained clinicians, as prescribed by a physician.
WARNING: RE-USE OF SINGLE-USE DEVICES CREATES A POTENTIAL RISK OF PATIENT OR USER. IT MAY LEAD TO CONTAMINATION AND/OR IMPAIRMENT OF FUNCTIONAL CAPABILITY. CONTAMINATION AND/OR LIMITED FUNCTIONALITY OF THE DEVICE MAY LEAD TO INJURY, ILLNESS OR DEATH OF THE PATIENT.
PRECAUTION: NOT INTENDED FOR USE WITH POWER INJECTORS.
DISPOSAL: Dispose of according to local guidelines and/or clinical protocol.

NOTICE TO THE USER: All serious incidents related to this product must be reported to B. Braun and the competent authority in the country where the product is being used.

INSTRUCTIONS FOR USE: Use Aseptic Technique.
• Visually inspect sterile barrier system for breaches and integrity prior to use. Do not use if packaging or product is damaged or contaminated or if protective caps are loose or missing. Use of damaged or contaminated device might lead to injury, illness or death of the patient.
• Replace set per institutional protocol.
• **To access CARESITE Luer Access Device:** Swab top of luer access device vigorously for 15 seconds with 70% Isopropyl Alcohol and allow to air dry prior to use. Use a luer connector without a needle. After medication administration, flush per institutional protocol, attaching completely between flushes.
1.Prime extension set to ensure air is expelled.
2.Firmly attach connectors to desired connections.
Unit label pictorial is for reference only.

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