



DENTAL SURGERY

COMPLETE PORTFOLIO FOR ORAL WOUND CLOSURE:
SUTURES AND BIOSURGICALS



COMPLETE CLOSURE TECHNOLOGIES PORTFOLIO FOR DENTAL SURGERY:

The needle-thread combinations listed in this catalog are just a selection of the several thousands possible combinations that exist in the B. Braun suture portfolio.

Therefore, if the desired combination is not listed here, please contact B. Braun Surgical, S.A.

Carretera de Terrassa, 121 | 08191 Rubí | Spain

Phone +34 93 586 6200 | Fax +34 93 588 1096 | www.bbraun.com

+++ Monosyn® Quick +++ Monosyn® +++ MonoPlus® +++ Novosyn® Quick +++
Novosyn® +++ Optilene® +++ Dafilon® +++ PremiCron® +++ Supramid +++ Silkam® +++

CONTENT

COMPLETE SUTURE RANGE

SUTURE MATERIAL NON-ABSORBABLE

SUTURE MATERIAL ABSORBABLE

B. BRAUN NEEDLE NOMENCLATURE BIOSURGICALS







CLOSURE TECHNOLOGIES SOLUTIONS

	4	
Supramid	6	
Dafilon®	8	
Optilene®	10	
Silkam®	12	
PremiCron®	14	
Monosyn® Quick	16	
Monosyn®	18	
MonoPlus®	20	
Novosyn® Quick	22	
Novosyn®	28	
	33	
Gelita	34	
Bone Wax	34	
Lyostypt®	35	
Cellistypt®	35	
	36	





COMPLETE SUTURE RANGE

NON-ABSORBABLE SUTURE MATERIAL

PSEUDO-MONOFILAMENT / MONOFILAMENT







Pseudo-monofilament Advantage	Monofilament Advantage	Monofilament Advantage
Supramid	Dafilon®	Optilene®
		
		
Polyamide 6/6.6, non-absorbable	Polyamide 6/6.6 (dyed blue or undyed), Polyamide 6.6 (dyed black), non-absorbable	Polypropylene and polyethylene, non-absorbable
Coating: none	Coating: none	Coating: none
USP 6/0 to 3 metric 0.7 to 6	USP 11/0 to 3 metric 0.1 to 6	USP 10/0 to 2 metric 0.2 to 5

BRAIDED





Braided	Braided
Silkam®	PremiCron®
	
	
Natural silk, non-absorbable	Polyester, non-absorbable
Coating: pure beeswax or refined paraffin wax	Coating: silicone
USP 8/0 to 6 metric 0.4 to 8	USP 6/0 to 5 metric 0.7 to 7

ABSORBABLE SUTURE MATERIAL


MONOFILAMENT

Monofilament Advantage	Monofilament Advantage	Monofilament Advantage
Monosyn® Quick	Monosyn®	MonoPlus®
		
		
Glyconate, short-term absorbable	Glyconate, mid-term absorbable	Polydioxanone, long-term absorbable
Coating: none	Coating: none	Coating: none
USP 6/0 to 1 metric 0.7 to 4	USP 7/0 to 2 metric 0.5 to 5	USP 7/0 to 2 metric 0.5 to 5


BRAIDED

Braided	Braided
Novosyn® Quick	Novosyn®
	
	
Polyglactin 910, short-term absorbable	Polyglactin 910, mid-term absorbable
Coating: polyglactin 370 and calcium stearate	Coating: polyglactin 370 and calcium stearate
USP 6/0 to 2 metric 0.7 to 5	USP 8/0 to 2 metric 0.4 to 5

Monofilament Advantage




Braided



Supramid









NON-ABSORBABLE SUTURE



DESCRIPTION

Supramid is a non-absorbable synthetic suture made of polyamide, available in two different structures depending on suture diameter:

- Monofilament (USP 6/0-5/0)
- Pseudo-monofilament (USP 4/0-2)









Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
3/8 circle reserve cutting needle									
	45 black			G0712078	G0712086	G0712094			
	45 black			G0712124	G0712132	G0712140			
	45 75 black			G0712191	G0712205 G0714208	G0712213 G0714216	G0712221 C0714224		
	45 75 black				G0712345	G0712353 G0714356	G0712361 G0714364	G0714372	
1/2 circle reverse cutting needle									
	45 black			G0712035	G0712043	G0712051			
	45 black 75 white			G0712248	G0712256	G0712264 C0692719			
	75 black						G0712426	G0712434	
	75 white						C0692689		

C Box of 36 units

B Box of 24 units

G Box of 12 units

M Box of 6 units










Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
5/8 circle conventional cutting needle									
	45 black				C0712884	C0712892			
	45 black					C0712930			
3/8 circle reverse cutting needle with precision point									
 	45 black		C0716073	C0716081					
1/2 circle reverse cutting needle with precision point									
 	45 black			C0716189					
 	45 black				C0716219				










	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
Pre-cut sutures									
5 threads of 45 cm	3 x 45 black			C0216445	C0216453				



DESCRIPTION


Dafilon® is a non-absorbable synthetic monofilament suture made of polyamide polymers 6/6.6 (dyed blue or undyed) or polyamide 6.6 (dyed black).

Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
3/8 circle reverse cutting needle									
	45 blue		G0932060	G0932078	G0932086				
	45 blue		G0932116	G0932124	G0932132				
	45 75 blue			G0932191	G0932205 G0935204	G0932213 G0935212	G0932221		
	45 75 blue				G0932345 G0935344	G0932353 G0935352	G0932361 G0935360		
	45 75 blue					C0932469 G0935468	G0935476		
3/8 circle trapezoidal micro engraved with precision point									
	45 blue			C0933908					
	45 blue			G0933906	G0933903				
	45 75 blue			G0933907	C0933909 C0933904				
	45 75 blue			G0933911	G0933910 C0933905	C0934008 C0933901			

Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
3/8 circle reverse cutting needle with precision point									
 DSMP 7 	45 blue		C0930018						
 DSMP 13 	45 blue		G0936081						
Needle length mm	Thread length cm color	USP (metric)							
		11/0 (0.1)	10/0 (0.2)	9/0 (0.3)	8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)
3/8 circle round bodied needle micro needle									
DRm 3 	5 15 black	G1118064	G1117041						
DRm 4 	10 black	G1118080	G1118099						
DRm 4s 	10 black			G1118137					
DRm 5 	15 black		G1118110	G1111434					
DRm 6 	15 black			G1118161					








Optilene®

NON-ABSORBABLE SUTURE



DESCRIPTION

Optilene® is a non-absorbable monofilament surgical suture made of polypropylene and polyethylene.






Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
3/8 circle reverse cutting needle									
 DS 12	45 blue		G3090206	G3090207					
 DS 16	45 75 blue		G3090211	G3090212 C3090512	G3090213 G3090513				
 DS 19	45 blue			G3090219	G3090220	G3090221			
 DS 24	45 75 blue				G3090434	G3090435	G3090436 G3090236		
 DS 30	45 75 blue					G3090246	C3090447 G3090247		
1/2 circle reverse cutting needle									
 HS 18	45 blue			C3090264					
 HS 26	75 blue					G3090269	G3090270		

C Box of 36 units

B Box of 24 units

G Box of 12 units










M Box of 6 units









Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
3/8 circle trapezoidal micro engraved with precision point									
DGMP 11 	45 blue		G3095588						
DGMP 13 	45 blue		G3095589	G3095518	G3095910				
DGMP 16 	45 blue			G3095549	C3095499				
DGMP 19 	45 blue			G3095550		C3095591			
3/8 circle reverse cutting needle with precision point									
DSMP 11 	45 blue	G3090606		G3090608					



DESCRIPTION

Silkam® is a non-absorbable, braided and coated suture, made of natural silk filaments.


Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
1/2 circle round bodied needle									
 HR 13	45 black			C0760072					
 HR 17	75 black				G0760137	G0760145			
 HR 22	75 black			C0760285		G0760307	G0760315		
 HR 26	75 black					G0760412	G0760420		
3/8 circle reverse cutting needle									
 DS 12	45 black		G0762067	G0762075	G0762083	G0762091			
 DS 16	45 75 black		G0762113	G0762121	G0762130	G0764140	G0764159		
 DS 19	45 75 black		C0765171	G0762199	G0762202	G0762210 G0765210	G0762229		
 DS 24	75 black				G0762342	G0762350	G0762369		
 DS 30	75 black						G0762474		G0762490

Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
1/2 circle reverse cutting needle									
<div>HS 15</div> <div></div>	45 75 black			C0765490	G0762040	G0762059			
<div>HS 21</div> <div></div>	45 black			G0762245	G0762253	G0762261			
<div>HS 23</div> <div></div>	45 black					G0762300	G0762326		
<div>HS 26</div> <div></div>	75 black					G0762415	G0762423		
3/8 circle reverse cutting needle with precision point									
<div>DSMP 7</div> <div><div></div></div>	45 black	G0766038							
<div>DSMP 11</div> <div><div></div></div>	45 black		G0766070	G0766089					

	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
Pre-cut sutures									
3 threads of 45 cm	3 x 45 black			C0266248	C0266256				
5 threads of 45 cm	5 x 45 black				C0266450				
10 threads of 45 cm	10 x 45 black			C0266647	C0266655	C0266663			
Ligature reels									
1 thread of 250 cm	1 x 250 black			G0260410	G0260428	G0260436			















PremiCron®

NON-ABSORBABLE SUTURE



DESCRIPTION

PremiCron® is a non-absorbable braided and coated synthetic multifilament composed of polyester fibres.




Needle length mm	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (3)	2 (5)
1/2 circle round bodied needle									
 	45 green		G0026029						
 	75 green			G0026004					
1/2 circle round bodied needle with short cutting point									
 	90 green				C2026305				
3/8 circle reverse cutting needle									
 	45 green		C0026172	G0026173					
	45 green		C0026176	G0026177					
 	45 green		C0026551	G0026552	G0026553	C0026554			
 	45 75 green			C0026514	C0026515	C0026190			
	75 green						G0026532	C0026533	

C Box of 36 units

B Box of 24 units

G Box of 12 units

M Box of 6 units

Needle length mm	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (3)	2 (5)
1/2 circle reverse cutting needle									
 HS 15	45 green			C0026197	C0026198				
 HS 21	45 green			C0026199	G0026200				
 HS 26	75 green				C0026420	C0026421			










Monosyn® Quick

ABSORBABLE SUTURE



DESCRIPTION

Monosyn® Quick is an undyed short-term absorbable synthetic monofilament suture made of glyconate for the approximation of soft tissue with a fast healing process.














Needle length mm	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (3)	2 (5)
3/8 circle reverse cutting needle									
 DS 12	45 70 undyed	C0025134	C0025135 C0025323	C0025136 C0025304	C0025169				
 DS 16	45 70 undyed		C0025125 C0025331	C0025126 C0025332	C0025333				
 DS 19	45 70 undyed	C0025005	C0025001 C0025006	G0025002 C0025007	C0025003 C0025008	C0025009			
 DS 24	45 70 undyed			C0025011	C0025351 C0025012	C0025352 C0025013	C0025342		
 DS 30	70 undyed				G0025014	G0025015	C0025107		
1/2 circle reverse cutting needle									
 HS 15	45 undyed		C0025143						
 HS 18	70 undyed				C0025379				
 HS 23	45 undyed			C0025388	C0025389				
 HS 37	70 undyed						C0025174		

C Box of 36 units

B Box of 24 units

G Box of 12 units















M Box of 6 units

Needle length mm	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (3)	2 (5)
1/2 circle reverse cutting needle									
	90 undyed					C0025265	C0025266		
	90 undyed						C0025256	C0025257	
3/8 circle trapezoidal micro engraved with precision point									
	45 undyed	G0025807	C0025805						
	45 undyed		G0025806	C0025802					
	45 undyed			G0025803					
	45 70 undyed		C0025808	G0025804	C0025801				
	45 70 undyed				C0025809 C0025138				
3/8 circle reverse cutting needle with precision point									
	70 undyed		C0025112						
	45 70 undyed	G0025023 C0025111	C0025113	C0025337					
	45 70 undyed	C0025027	G0025028	C0025130	G0025030				
	45 70 undyed			C0025036	C0025033	G0025038			
	45 70 undyed			C0025039 C0025042		C0025044			
1/2 circle reverse cutting needle with precision point									
	45 undyed			C0025101					













DESCRIPTION

Monosyn® is a mid-term absorbable synthetic monofilament suture made of glyconate, available in undyed color for the closure of superficial wounds and in violet color for the rest of indications. Monosyn® offers high knot pull tensile strength.

Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
3/8 circle reverse cutting needle									
 	45		G0023102	G0023103					
	70		G0023412	G0023413	C0023414				
	70		G2023412	G2023413	C2023414				
	undyed								
 	45			G0022419					
	45			G2022419					
	violet								
 	45		C0023562	C0023563	G0023303				
	45				G2023303				
	70				G0023424	G0023425			
	70				G2023424	G2023425			
	undyed								
 	70				G0022204	G0022205			
	70				G2022204	G2022205			
	violet								
 	45			C0023503	G0023204	G0023205			
	45				G2023204	G2023205			
	70			G0023403	G0023404	G0023405	C0023406		
	70			G2023403	G2023404	G2023405			
	undyed								
 	70						G0022216		
	violet								
	45					C0023115	C0023116		
 	70				G0023214	G0023215	G0023216		
	70				G2023214	G2023215			
	undyed								

Italic codes: RCP packaging

Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
1/2 circle reverse cutting needle									
 HS 18	70 undyed			C0023263	C0023264				
 HS 23	70 70 undyed				C0023224	G0023225 G2023225	C0023226		
 HS 26	70 violet 70 undyed						G0022237 G0023240		
 HS 26s	70 violet						G0022242	G0022243	
3/8 circle trapezoidal micro engraved with precision point									
 DGMP 11	45 undyed		C0023709	C0023705					
 DGMP 13	45 undyed		C0023710	G0023706	G0023702				
 DGMP 16	45 70 undyed			C0023707	G0023703	C0023711			
 DGMP 19	45 45 70 70 undyed			C0023708	G0023704 G2023704	G0023701 G2023701			
 DGMP 24	70 70 undyed					C0023717 C2023717			
3/8 circle reverse cutting needle with precision point									
 DSMP 19	70 70 undyed				G0023627 G2023627				

Italic codes: RCP packaging













MonoPlus®

ABSORBABLE SUTURE



DESCRIPTION

MonoPlus® is a violet long-term absorbable synthetic monofilament suture made of polydioxanone. MonoPlus® offers good handling combined with great strength. Due to its slow degradation profile, MonoPlus® is indicated for use in general soft tissue approximation, specially in cases where an extended wound support of more than 4 weeks is desirable.

Needle length mm	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 HR 13ss 	70 violet				C0024205				
 HR 17 	70 violet		G0024013	G0024014					
 HR 22 	70 violet			G0024134	G0024135				
 HR 22s 	70 violet						C0024147		
 HR 26 	70 violet			G0024024	G0024025	G0024026			
1/2 circle reverse cutting needle									
 HS 26s 	70 violet					G0024186			

C Box of 36 units

B Box of 24 units

G Box of 12 units

M Box of 6 units

















Novosyn® Quick

ABSORBABLE SUTURE



DESCRIPTION

Novosyn® Quick is a polyglactin 910 braided suture with fast degradation profile. Novosyn® Quick is intended for use in soft tissue approximation when only short-term wound support is required and when the rapid absorption of the suture would be beneficial.











Needle length mm	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 	70 undyed		C3046002						
 	70 undyed		C3046012	C3046013	C3046014				
 	45 70 undyed			C3046029	C3046024 C3046030	C3046031			
 	70 undyed					C3046131			
 	45 70 undyed			C3046040	C3046051 C3046041	C3046042	C3046043		
 	70 undyed					C3046090	C3046091	C3046092	
 	70 undyed				C3046046	C3046047	C3046048	C3046049	
 	70 undyed							C3046053	

C Box of 36 units

B Box of 24 units










G Box of 12 units

M Box of 6 units

Needle length mm	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
	70 undyed						C3046006		
	45 90 undyed					C3046590 C3046595	C3046596	C3046597	
	70 90 undyed				C3046554	C3046467	C3046556	C3046557	C3046559 C3046558
	90 undyed					C3046170			
	90 undyed						C3046506		
	90 undyed					C3046180	C3046661	C3046663	
5/8 circle round bodied needle									
	90 undyed						C3046083		
	90 undyed						C3046078		
3/8 circle round bodied needle with short cutting point									
	70 undyed		C3046311						
	70 undyed			C3046320	C3046321				

Novosyn[®] Quick

ABSORBABLE SUTURE










Needle length mm	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle with short cutting point									
	90 undyed					C3046651	C3046652	C3046653	
	90 undyed					C3046664	C3046662	C3046365	
	90 undyed							C3046375	
1/2 circle round bodied needle with blunt point									
	90 undyed						B3046605		
3/8 circle reverse cutting needle									
	70 undyed		C3046202	C3046203					
	70 undyed		C3046212	C3046213	C3046214				
	45 70 undyed		C3046219 C3046419	C3046220 C3046420	C3046221 C3046421	C3046422			
	45 70 undyed			C3046234	C3046535 C3046235	C3046536 C3046236			
	90 undyed				C3046246	C3046247	C3046248		

C Box of 36 units

B Box of 24 units






G Box of 12 units

M Box of 6 units

Needle length mm	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle reverse cutting needle									
 HS 15	45 undyed			C3046204					
 HS 18	45 70 undyed			C3046226	C3046227				
 HS 23	45 undyed			C3046229	C3046230	C3046231			
 HS 37s	90 undyed					C3046240	C3046241	C3046242	
 HS 40s	90 undyed							C3046587	
 HS 48	90 undyed						C3046525		
Straight reverse cutting needle									
 GS 13	45 undyed			C3046208					
 GS 51	70 undyed			C3046271	C3046272	C3046273	C3046274		
 GS 60	70 undyed					C3046410			

Novosyn® Quick

ABSORBABLE SUTURE

Needle length mm	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
3/8 circle reverse cutting needle with precision point									
 DSMP 11	45 70 undyed	C3046607 C3046603	C3046610						
 DSMP 13	45 undyed		C3046619	C3046620					
 DSMP 16	45 70 undyed		C3046611	C3046625 C3046685					
 DSMP 19	45 70 undyed			C3046629 C3046529	C3046630 C3046530	C3046631 C3046531			
 DSMP 24	45 70 undyed				C3046154 C3046657	C3046155 C3046658			


C Box of 36 units

B Box of 24 units

G Box of 12 units



















M Box of 6 units

	Thread length cm color	USP (metric)							
		5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	3 (6)
Pre-cut sutures									
3 threads of 45 cm	3 x 45 undyed			C3037625	C3037626	C3037627	C3037628		
6 threads of 45 cm	6 x 45 undyed		C3037004	C3037005	C3037006	B3037007	B3037008		
2 threads of 70 cm	2 x 70 undyed				C3037026	C3037027	C3037028		
5 threads of 70 cm	5 x 70 undyed		C3037444		C3037446				



DESCRIPTION

Novosyn® is a mid-term absorbable braided and coated synthetic suture made of polyglactin 910, that degrades by hydrolysis.























Needle length mm	Thread length cm color	USP (metric)							
		8/0 (0.4)	7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)
3/8 circle round bodied needle									
 DR 12 	60 70 violet			C0068434	C0068435	C0068105 C0068436			
 2xDR 12 	45 undyed			C0069725					
 DR 15 	45 70 violet					C0068137 C0068107			
 DR 18 	45 70 violet				C0068155	C0068138 C0068150	C0068151		
 DR 22 	70 violet					C0068159	C0068160	C0068161	
 DR 22b 	70 violet					C0068108	C0068117		
 DR 26 	70 120 violet						C0068164 C0068179	C0068165	
 DR 32 	70 violet						C0068170	C0068171	C0068172
Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
 DR 32 	70 90 violet					C0068566	C0068567	C0068173 C0068568	C0068168



















C Box of 36 units

















B Box of 24 units







G Box of 12 units

M Box of 6 units

Needle length mm		Thread length cm color	USP (metric)							
			7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
1/2 circle round bodied needle										
	HR 6 	30 45 violet	C0068149 C0068020							
	HR 10 	70 violet 45 undyed		C0068002 C0069802	C0068003					
	HR 13 	45 70 90 violet 45 70 undyed		C0068449 C0068006 C0069405	C0068450 G0068007 C0069407 C0069007	C0068451 C0068008 C0068452 C0069408 C0069008	C0068009			
	HR 13  take off	4 x 45 8 x 45 violet			C0088718 M0088190	C0088719 M0088191				
	HR 13ss 	45 70 violet					C0088320 C0068120	C0088321 C0068110	C0088322 C0068111	
	HR 17f 	70 violet				C0068888				
	HR 17 	45 70 90 violet 70 undyed		C0068011	G0068012 C0069012	C0088340 G0068013 C0068853 C0069013	G0068014 C0068854 C0069014	C0068015 C0068855 C0069015	C0068016	
	HR 17b 	70 violet			C0068144	C0068147	C0068152			
	HR 17  take off	4 x 45 8 x 45 4 x 70 8 x 70 violet 4 x 45 undyed				C0088713 M0088185 C0088613 M0088255 C0089713	C0088714 M0088187 C0088614			
	HR 17b  take off	8 x 45 violet						M0068153		
	2xHR 17 	70 90 violet			C0068869	C0068870	C0068904			

Needle length mm	Thread length cm color	USP (metric)							
		7/0 (0.5)	6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)
1/2 circle round bodied needle									
 2xHR 17b 	70 violet			C0068163	C0068176				
 HR 17s 	45 70 violet				C0068113	C0088335 C0068114	C0068115	C0068116	C0068118
 HR 20 	70 90 violet				C0068765	C0068766 C0068763	C0068767		
 HR 22 	20 45 70 90 120 violet 70 90 undyed			C0068028	G0068029 C0068438	C0068069 C0088325 G0068030 C0068877 C0068078	C0068004 C0068031 C0068881 C0068880	C0068032	
 HR 22b 	70 120 violet				C0068181	C0068186 C0068365	C0068189		
 HR 22  take off	4 x 45 8 x 45 4 x 70 8 x 70 violet				C0088729 M0088177 C0088629 M0088248	C0088730 M0088178 C0088630 M0088249	C0088731 M0088180 B0088631 M0088250		
 HR 22b  take off	4 x 45 8 x 45 violet				M0068182	C0088732 M0068188			
 2xHR 22 	90 violet					C0068874	C0068875		
 2xHR 22b 	70 90 violet				C0068190 C0068191	C0068192			

Needle length mm	Thread length cm color	USP (metric)								
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)	
1/2 circle round bodied needle										
 	70 violet 70 undyed			C0068932	C0068934 C0069934	C0068036	C0068037	C0068038		
 	70 violet						C0068937			
 	70 violet							C0068033		
 	20		C0068039	G0068040	C0068079	C0068018	C0068043	C0068044	C0068034	
	45				C0068109					
	70									
	4 x 70				C0068041	C0068042				
	90				C0088647					
	120				C0068540	C0068541				C0068542
 	violet				C0068088	C0068067				
	70			C0069040	C0069041	C0069042	C0069043			
	90				C0069541	C0069540	C0069543			
	undyed									
	 	70 violet			C0068196	C0068198	C0068209			
		4 x 45			C0088740	C0088741	C0088742			
8 x 45				M0088170	M0088171	M0088172				
 	4 x 70			C0088640	C0088641	C0088642				
	8 x 70			M0088237	M0088239	M0088241				
	violet									
	4 x 45				C0089741					
	8 x 45				M0089116	M0089117				
	undyed									
 	8 x 45 8 x 70 violet				M0068200 M0068201					

Needle length mm	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
1/2 circle round bodied needle									
 2xHR 26 	70 90 violet				C0068901	C0068902	C0068907		
 2xHR 26b 	70 90 violet				C0068210 C0068216	C0068217			
 HR 26s 	45 violet					C0068019			

	Thread length cm color	USP (metric)							
		6/0 (0.7)	5/0 (1)	4/0 (1.5)	3/0 (2)	2/0 (3)	0 (3.5)	1 (4)	2 (5)
Pre-cut sutures									
3 threads of 45 cm	3 x 45 violet			G0058624	G0058625				

B. BRAUN NEEDLE NOMENCLATURE

TO IDENTIFY THE DIFFERENT TYPES OF NEEDLES, B. BRAUN USES THE FOLLOWING NOMENCLATURE:


A	B	C	12	D
Needle curvature	Body type	Point type	Needle length	Remark
S 1/8 circle	R round body	T trocar point	This refers to the distance between the needle point and the needle end, measured along the length of the needle in mm.	ff extra thin bodied needle
V 1/4 circle	S cutting needle	N blunt point		f thin bodied needle
D 3/8 circle	L lancet needle	S sternum point		s heavy bodied needle
H 1/2 circle	G trapezoidal body	C short cutting point		ss extra heavy bodied needle
F 5/8 circle		M micro-point		v break-off needle
G straight		MP precision point		b black needle
P progressive curved		m micro needle		
J hook		mV micro vascular		
SK ski needle		B protection point		

Example:

H	R	C	26	
----------	----------	----------	-----------	--

HRC 26 → 1/2 circle round bodied needle with short cutting point, 26 mm long.

Needle angle of curvature




- S 1/8 circle
- V 1/4 circle
- D 3/8 circle
- H 1/2 circle
- F 5/8 circle

"Take off" removable needle

In order to remove the needle from the suture without the aid of scissors, one must pull the suture taut and secure and grasp the needle in the needleholder. A straight tug of the needleholder will then release the needle from the suture.

Others



BIOSURGICALS



Size	Content/box	Code
1 x 1 x 1 cm	10 tampons	2070014
1 x 1 x 1 cm	50 tampons	2070103
1 x 1.5 x 1.5 cm	50 tampons	2070154
8 x 2 x 1 cm	10 tampons	2070707
8 x 5 x 1 cm	10 tampons	2070600

Gelita

- Moldable and compressible (1)
- High absorption capacity: water uptake ≥ 35 fold its own weight (2)
- Absorbed within 4-6 weeks (3, 4)
- Not tissue-irritating (5)
- Can be used as carrier for thrombin (6, 7)



Bone Wax

- Mechanical haemostasis in bones (8-10)
- Haemostasis achieved by tamponade effect (8-10)
- Achieves mechanical occlusion of intraosseous vessels (8-10)
- Soft, easy to shape and apply "only hand-warm" (11-13)

Format	Content/box	Code
Plate	24 pieces	1029754



Lyostypt®

- Local hemostatic agent for minor bleeding
- Absorbable collagen fleece
- Absorbable within 3-4 weeks

Size	Content/box	Code
3 x 5 cm	12 pieces	1069128
5 x 8 cm	6 pieces	1069152
5 x 8 cm	12 pieces	1069020
10 x 12 cm	4 pieces	1069209
5 x 30 cm	4 pieces	1069306
10 x 12 cm	8 pieces	1069039



Cellistypt®

- Bacteriostatic properties (11, 14, 15)
- Absorbed in 1-2 weeks (16)
- Easy to use (16)
- Biocompatible (17)

Cellistypt®

Size	Content/box	Code
5 x 1.25 cm	15 pieces	2080501
5 x 7 cm	15 pieces	2080508
7 x 10 cm	15 pieces	2080511
5 x 35 cm	10 pieces	2080536
10 x 20 cm	10 pieces	2080541
1.5 x 1.5 cm	40 pieces	2080515

Cellistypt® D-K

Size	Content/box	Code
2.5 x 2.5 cm	15 pieces	2081203
2.5 x 9 cm	15 pieces	2081209
5 x 7.5 cm	10 pieces	2081275
7 x 10 cm	10 pieces	2081210
14 x 20 cm	10 pieces	2081240

Cellistypt® F

Size	Content/box	Code
2.5 x 5 cm	10 pieces	2082025
5 x 7.5 cm	10 pieces	2082075
5 x 10 cm	10 pieces	2082005
10 x 10 cm	10 pieces	2082010
10 x 20 cm	10 pieces	2082020

CLOSURE TECHNOLOGIES SOLUTIONS

By 1906, more than a hundred methods for sterilizing catgut existed. The breakthrough came in 1908, when a pioneering company created a technique able to assure the supply of consistently sterile sutures to hospitals and other medical facilities. This company was B. Braun!

Since then, B. Braun is committed to innovation and excellence in product-driven research, development and production. Furthermore, in Closure Technologies (CT) we remain loyal to our roots by working closely together with healthcare professionals, making them part of our principal commitment. In CT we are *Closer to You*.

COMPLETE PRODUCT RANGE

B. Braun offers a complete range of absorbable sutures, non-absorbable sutures, customized surgical procedures sutures sets (Pro-Set) and special suture sets, as well as a full range of meshes, biosurgicals and other products for wound closure.

HIGH QUALITY

B. Braun has achieved its leading position in the field of surgery thanks to its high quality products. All B. Braun manufacturing plants have a certified quality management system for assuring the quality of raw materials, semi-finished and finished products, and guarantees the traceability of each product manufactured.

INNOVATIVE

It is our belief that only innovative products, which increase the customer benefit, can ensure a long term success on the market. Our innovative products, developed by our own R&D department, are also the result of the shared knowledge between competent employees as well as clinical experts and users, true to our motto: Sharing Expertise.

AESCULAP ACADEMY

The Aesculap Academy is one of the leading medical forums for everyone who is professionally, passionately and ambitiously committed to people's health. To those medical professionals we offer top quality knowledge transfer based on globally recognized quality criteria using innovative methods and technologies.



REFERENCES

- (1) Determination of the compression and mouldable properties of Gelita Tampon according test RDR-DID-AJR-16191.
- (2) Determination of water absorption for Gelita Tampon according test RDR-DID-MLB-AJR-16192.
- (3) Broekema FI, Van Leeuwen MB, Van Minnen B, Bos RR. IN vivo degradation of polyurethane foam with 55 wt % polyethylene glycol. *J Biomed Mater Res A*. 2015 Nov;103(11):3666-75.
- (4) Mani M, Ebenezer V, Balakrishnan R. Impact of Hemostatic Agents in Oral Surgery. *Biomed Pharmacol J*. 2014 Feb 06;7(1):215-19.
- (5) Based on biocompatibility tests performed according to the requirements of EN ISO 10993-1:2009/AC:2010.
- (6) Yao HH, Hong MK, Drummond KJ. Haemostasis in neurosurgery: what is the evidence for gelatin-thrombin matrix sealant? *J Clin Neurosci*. 2013 Mar;20(3):349-56.
- (7) Masci E, Santoleri L, Belloni F, Bottero L, Stefanini P, Faillace G, Bertani G, Montinaro C, Mancini L, Longoni M. Topical hemostatic agents in surgical practice. *Transfus Apher Sci*. 2011 Dec;45(3):305-11.
- (8) Stehrer R, Hunger S, Schotten KJ, Parsaei B, Malek M, Jacob M, Geiselseder G, Meier J. Reduction of transfusion requirements in pediatric craniostylosis surgery by a new local hemostatic agent. *J Craniomaxillofac Surg*. 2016;44(9):1246-51.
- (9) Koc K, Ilik MK. Surgical management of an osteoblastoma involving the entire C2 vertebra and a review of literature. *European Spine Journal*. 2016;25 Supplement 1:220-3.
- (10) Elmorsy A, Abdel-Wahab A, Sabry M, Wahby E. Bone wax versus oxidized regenerated cellulose in preventing sternal intramedullary bleeding in cardiac surgery. *Journal of the Egyptian Society of Cardio-Thoracic Surgery*. 2016;24(4):302-7.
- (11) Schonauer C, Tessitore E, Barbagallo G, Albanese V, Moraci A. The use of local agents: bone wax, gelatin, collagen, oxidized cellulose. *Eur Spine J*. 2004;13(Suppl. 1):89-96.
- (12) Ogle OE. Perioperative hemorrhage. In: Dym H, Ogle OE. *Atlas of Minor Oral Surgery*. Philadelphia, Pa: WB Saunders. 2000:62-3.
- (13) Santhosh Kumar MP. Local hemostatic agents in the management of bleeding in oral surgery. *Asian Journal of Pharmaceutical and Clinical Research*. 2016;9:3.
- (14) Spangler D, Rothenburger S, Nguyen K, Jampani H, Weiss S, Bhende S. In vitro antimicrobial activity of oxidized regenerated cellulose against antibiotic-resistant microorganisms. *Surg Infect (Larchmt)*. 2003;4(3):225-62.
- (15) Wagner WR, Pachence JM, Ristich J, Johnson PC. Comparative in vitro analysis of topical hemostatic agents. *J Surg Res*. 1996;66:100-8.
- (16) See Cellistyp IFU.
- (17) Data on file, report AP-SY-1501.

Distributed by:

B. Braun Surgical, S.A. | Carretera de Terrassa, 121 | 08191 Rubí | Spain

Phone +34 93 586 6200 | Fax +34 93 588 1096 | www.bbraun.com

AESCULAP[®] – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany

Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

Manufacturer acc. to MDD 93/42/EEC

B. Braun Surgical, S.A. | Carretera de Terrassa, 121 | 08191 Rubí | Spain

The main product trademark "Aesculap" and the product trademarks "Cellistyp", "Dafilon", "Lyostyp", "MonoPlus", "Monosyn", "Novosyn", "Optilene", "PremiCron" and "Silkam" are registered trademarks of B. Braun Melsungen AG and its subsidiaries.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.